

EYECAP™
Eye Irrigation Shield

E

Y

E

C

A

P

D O N T W A I T

T O I R R I G A T E

- Instant Eye Irrigation
- Simple & Easy to Use
- Avoids Blinking Reflex
- Fits Any Standard Saline Bottle
- Child & Adult Friendly



EYECAP™

Eye Irrigation Shield

EyeCap™ is a twist-on eye irrigation shield that allows fast and easy eye irrigation.

Use EyeCap™ for immediate eye irrigation at triage, by EMS or for occupational splash injuries. The EyeCap™ can serve as a bridge to other irrigation methods that take longer to set up and are difficult for most patients to tolerate.

How it Works

The wide contoured rim of the EyeCap™ seals easily over the orbit. The nozzle creates a gentle, low pressure spray that deflects off the sidewalls and trickles down to create a “pool” of irrigation fluid. Now patients can easily open eyelids as if they were swimming underwater. An overflow port drains run-off and allows continuous flow of high volumes of fluid.

Benefits

- Don't wait to irrigate
- Fits any standard saline bottle
- Simple to use
- Avoids triggering blinking reflex
- Instantly begin dilution treatment
- Eye remains open for better surface exposure
- High volume continuous dilution
- Comfortable for patient
- Less staff required

Ordering Information

Code Number	Description	Qty.
7446802	Eye Irrigation Shield	25 / box



Please refer to the product label and / or package insert for full instructions on the safe use of these products.

EyeCap™ is a trademark of Splash Medical Devices, LLC. Patented / Patent Pending.


Distributed by

BSN medical Inc. 5825 Carnegie Blvd. Charlotte, NC 28209-4633
Tel. (+1) 704 554 9933 Fax (+1) 704 358 4558 www.bsnmedical.com
To order toll-free: BSN medical (+1) 800 552 1157



60636 R1

© 2011 BSN medical Inc.



EYECAP

REV 05/11