

EYECAP™

EYE IRRIGATION ADAPTER



- Fits on Saline Bottles
- Instant Treatment
- Simple & Fast
- Non-Invasive
- Portable
- Pooled Dilution
- Kid Friendly

EYECAP™

EYE IRRIGATION ADAPTER

Quick & Non-Invasive

Zero Impact Spray Nozzle

Patient Opens Eye in "Swimming Pool"

Overflow Port Allows Continuous Irrigation



Multiple Uses:

- Splash Injuries
- Foreign Bodies
- Allergens
- Irritants
- Decon

Zero Impact Nozzle:

- Stream Aims Outward
- Deflects Off Side Wall
- Trickles Down
- Gently Pools Over Orbit
- Avoids Blink Reflex

Swimming Pool Effect:

- Fluid Retained By Seal
- Large Volume Dilution
- Wide Comfort Rim
- Patient Opens Eye In Underwater Bath

Continuous Irrigation:

- Fluid Exits Sideport
- Ongoing Dilution Effect While Seal Maintained.

DON'T WAIT TO IRRIGATE! Use the EyeCap for quick and immediate irrigation of the eye. It is Simple and Easy. Use the EyeCap as a bridge to irrigation lens that may take longer to set up and can be difficult to tolerate. Use for definitive treatment of common minor injuries. The EyeCap fits on compatible saline bottles. It is compact and portable for use on scene, on the job, en route care, at triage and in an exam room. It is non-invasive and kid friendly. Begin irrigation in seconds. Stop the Sizzle! Flush out irritating or abrasive foreign bodies.

Just twist on any compatible saline bottle and squeeze. The zero impact spray nozzle hits the sidewall. Fluid then trickles down and pools up over the orbit to comfortably bathe the eye surface and soak under the lids. Refer to website or packaging for additional info. Non-Sterile. For single use only. Do not use if a ruptured globe is known or suspected.



Seal EyeCap Over Orbit.
Squeeze Irrigation Bottle.
Fluid Pools Over Eye.
Start Irrigation Instantly.

Decontaminate Eye
External Surfaces by
Releasing Initial Fluid.
Next Keep Fluid Sealed.

Assist Opening Lids Initially.
Fluid Penetrates Under Lids
For Quick Dilution and Relief.
Relief Enhances Cooperation.

Large Volume Irrigation Pool.
Contain Fluid Within Adapter
Overflow Exits via Side Port.
Continue With More Fluid.

EyeCap Can Be Definitive Treatment.
Simple & Instant First Line Treatment.
Stop the Sizzle! Don't Wait to Irrigate!
IF Still Indicated, Insert Irrigation Lens.