

SPLASHCAP^M

Wound Irrigation Shield



- Faster Than Syringes
- Universal Fitting
- Built-in Shield
- Kid Friendly
- Economical



SPLASHCAP[™]

Wound Irrigation Shield

SPLASHCAP[™] is a twist-on wound irrigation shield that allows a clinician to deliver hundreds of milliliters of irrigation fluid directly to a wound in seconds.

Features & Benefits

- Faster Wound Irrigation: Improved productivity and faster turnaround times.
- Safer Built-in Shield: Reduces risk of needle sticks and splattering when using bottle caps without shields.
- Proven Pressure (PSI): 4 15 PSI.
- Universal Fitting: Fits all major standard saline bottles.
- Better Wound Care: Standardize wound irrigation procedures in your department.
- High Volume: Deliver over 250 ml of irrigation fluid in seconds.
- Kid Friendly: No scary syringes.



SPLASHCAP™ advantage over a syringe shield

- Time savings
- Easier to use
- · No repetitive drawing up the solution and flushing the wound
- · Attach and disconnect the shield to syringe each time
- A 20cc syringe = 25 times to deliver 500cc SPLASHCAP™ can deliver that volume in seconds

SPLASHCAP™ advantage over punching a hole in saline cap

- Eliminates needle-sticks
- Eliminates splatter from wound site
- Less mess

SPLASHCAP™ advantage over shields that require their saline bottle

- Reduced inventory
- Less expensive
- · Less storage space required



Code Number	Description	Qty.
7446800	SPLASHCAP™ Wound Irrigation Shield	25/case
7446801	SPLASHCAP™ Wound Irrigation Shield	100/case









SPLASHCAP™ is a trademark of Splash Medical Devices, LLC.



BSN medical Inc. 5825 Carnegie Blvd. Charlotte, NC 28209-4633 Tel. (+1) 704 554 9933 Fax (+1) 704 358 4558 www.bsnmedical.com To order toll-free: BSN medical (+1) 800 552 1157

60535 RN © 2010 BSN medical Inc. Printed in U.S.A. REV 01/10