


BABYSUCTION™

PEDIATRIC NASAL SUCTION



- Pediatric
- Nasal Adapter
- Blunt Tip
- One-Handed

BABYSUCTION™

BLUNT TIPPED PEDIATRIC NASAL SUCTION



The BabySuction is a pediatric suction device with a blunt nasal tip. The tip is designed to more easily create a seal with the outer surface of the nose when used with vacuum powered suctioning and to reduce trauma to the delicate internal surfaces of the nose. Upgrade to the TwinSuction for an even more versatile device. Consider using it with nasal irrigation fluid for an improved nasal lavage technique. Non-sterile. For single patient use only. For single use only & professional use only.

Nasal Suction

- Pediatric Suction
- Blunt Nasal Tip Adapter
- Mechanical Vacuum Device

Reduce Trauma Risks

Simpler Set Up

Improved Preparation