

EASIEAR™

SMOOTH WIRE CURETTE



- Smooth Round Wire
- No Abrasive Edges
- Gentle Contact
- Flexible Handle
- Rigid Tip
- Disposable
- Inexpensive

EASIEAR™

SMOOTH WIRE CURETTE



Use the EasiEar to remove wax and foreign bodies. The smooth round wire tip has no abrasive edges. The spring flexes to accommodate surface contours. The rigid angled tip allows you to effectively scoop out wax and foreign bodies from the ear. It can even be used to scoop out foreign bodies from the nose!

Not only is the EasiEar a great product, it **costs less than plastic cures**. Also, consult with your billing professional regarding the potential use of reimbursement codes when performing ear wax and foreign body removal.

COMPARE TIP SURFACES:

The EasiEar™ round wire curette has a smooth gentle tip surface. Compare this to some earlier injection molded plastic cures. Notice the abrasive flash at the edges of the parting lines, where the mold plates separated. That can't feel good scraping on the sensitive ear canal. You may not typically see this, but your patients are feeling it. Try our product and help us: Make it EasiEar for Everyone!

