


# TWINSUCTION™


## REVERSIBLE PEDIATRIC SUCTION



- Pediatric
- Reversible
- 2 in 1 Device
- Blunt Nasal Adapter
- Oropharyngeal Catheter
- Soft Flexible Tubing
- Simpler Stocking & Prepration

# TWINSUCTION™

## REVERSIBLE PEDIATRIC SUCTION



The TwinSuction is a reversible 2-in-1 pediatric suction device for use with mechanical suction. One end has a nasal tip to more easily create a seal with the outer surface of the nose when used with vacuum powered suctioning. The blunt shape is designed to reduce trauma to the delicate internal surfaces of the nose. The opposite end has a soft flexible catheter for deeper oropharyngeal suctioning. Now it is possible to stock a single pediatric device that performs both functions. Using the TwinSuction may also avoid the need for deep suctioning and its associated complications. Of course continue to use that technique when indicated.

Choose the end desired. Attach to a vacuum hose. If your needs change, just flip it!  
Non-sterile. For single patient use only. For single use only & professional use only.



### Reversible Suction

- 2 in 1
- Pediatric Suction
- Blunt Nasal Tip Adapter
- Soft Flexible Oropharyngeal

Reduce Trauma Risks

Simpler Set Up

Improved Preparation