

BABYSUCTION™


PEDIATRIC NASAL SUCTION



- Pediatric
- Nasal Adapter
- Blunt Tip
- One-Handed

BABYSUCTION™

BLUNT TIPPED PEDIATRIC NASAL SUCTION



The BabySuction is a pediatric suction device with a blunt nasal tip. The tip is designed to easily create a seal with the opening of the nose when vacuum powered suctioning is used. The device is also designed to reduce trauma to the delicate internal surfaces of the nose since the wide tip prevents it from going in too deeply. The rigid device can easily be manipulated with one hand. This leaves the other hand free for positioning and stabilizing the patient during the suctioning procedure.

Upgrade your equipment to the BabySuction for an easier to use device. Consider using it with nasal irrigation fluid in order to improve the nasal lavage technique. Non-sterile. For single patient use only. For single use only & professional use only.

Nasal Suction

- Pediatric Suction
- Blunt Nasal Tip Adapter
- Mechanical Vacuum Device

Reduce Trauma Risks

Simpler Set Up

One Handed Use

Improved Preparation